

Emergency Trolley Check List

On Top

Zoll defib with pads attached

On Drip Stand

ADULT Gum Bougee 15ch x60cm

Ambu-Bag with mask attached (inside bag)

Stethoscope

Drawer 1 AIRWAY

mac blade size 3 and 4 – both attached to laryngoscope handles
(check light can be activated on both laryngoscopes)

Macgills forceps

Oral airway size 2,3 & 4

Oral airway size 6 & 7

ET tube size 6, 7, 8 & 9

Stylet size medium and large

Catheter mount

ventilator filter

PEEP valve

Scalpel

Scissors

Et tube tie

Swabs

20ml syringe (prepared with Swab inside)

Aqualgel

Drawer 2 EMERGENCY DRUGS

X4 Adrenaline pre filled syringes 1:10000

x1 Amiodarone pre frilled syringe 300mg in 10ml

x1 Atropine pre filled syringe 1mg in 5ml

x1 10% Calcium chloride pre filled syringe

x5 Atropine 600mcg/1ml ampoules

x4 Naloxone 400mcg/1ml ampoules

x2 Chlorphenamine 10mg/1ml ampoules

x5 Adrenaline 1:1000 ampoules

X1 Dextrose 20% in 100mls

Air inlet

Giving set

5x Sodium Chlorine 0.9% 10ml plastic ampoules

Syringes 2ml, 5ml, 10ml, 20ml (5 of each)

Needles blue, green & white

Drawer 3 IV ACCESS

Tourniquet x2
x10 mediswabs
Venflons 2 of each, blue, pink, grey & orange
Veneguards x6
Needles x5 blue & green
Syringes x2 10ml, 2x 5ml
venflon dressing x3
x2 double lumen octopus
ABG syringe x4
x2 bandage
Micropore tape
Swabs

Drawer 4

500ml 0.9% Sodium Chloride x2
500ml Hartmans x2
Giving Set x4
Defib pads
Cardiac monitoring electrodes
Big Fish Ligature cutter
Flexible suction catheters size 12 + 14
Yanker for suction

Side of trolley

Full Oxygen cylinder with suction unit assembled, tubing and yanker attached
Sharps box

Defibrillator

The defib auto checks at midnight a green tick can be seen on window to the top right, If any faults occur a red cross will show at the top right.

If a fault occurs check all connections, check defib pads are attached and carry out a manual test, instructions attached.

If a green tick is not seen report to medical physics.

Oxygen Cylinder

This will empty in 8 minutes if used to deliver high flow O2 and suction, replace cylinder if used

Cleaning

Please clean the trolley as part of the daily checks, including the front & back of the trolley

Covid-19 Precautions

A trolley set up nearby should contain;

x5 surgical gown

x5 visors

x5 surgical hat

a selection of the different FFP3 masks available